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journal homepage: [www.elsevier.com/math](http://www.elsevier.com/math)

Review article

Effectiveness of proprioceptive exercises for ankle ligament injury in adults:  
A systematic literature and meta-analysis<sup>☆</sup>

K. Postle<sup>a,\*</sup>, D. Pak<sup>b</sup>, T.O. Smith<sup>a</sup>

<sup>a</sup> School of Allied Health Professions, Faculty of Medicine and Health Science, University of East Anglia, Queen's Building, Norwich NR4 7TJ, UK

<sup>b</sup> Norwich Medical School, Faculty of Medicine and Health Science, University of East Anglia, Norwich NR4 7TJ, UK

“ The results indicated that **there is no statistically significant difference in recurrent injury** between the addition of proprioceptive exercises during rehabilitation of patients following ankle ligament injury “

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
### ORIGINAL RESEARCH ARTICLE

#### Ankle Muscular Proprioceptive Signals' Relevance for Balance Control on Various Support Surfaces An Exploratory Study

#### ABSTRACT

Forestier N, Terrier R, Teasdale N: Ankle muscular proprioceptive signals' relevance for balance control on various support surfaces: an exploratory study. Am J Phys Med Rehabil 2014;00:00-00.

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## Impaired control of weight bearing ankle inversion in subjects with chronic ankle instability

R. Terrier<sup>a,b,\*</sup>, K. Rose-Dulcina<sup>b</sup>, B. Toschi<sup>b</sup>, N. Forestier<sup>b</sup>

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### ARTICLE INFO

Article history:  
Received 4 September 2013  
Accepted 9 January 2014

Keywords:  
Ankle  
Inversion  
Sprain  
Chronic instability  
Evertor muscle weakness  
Weight bearing  
Angular velocity

### ABSTRACT

**Background:** Previous studies have proposed that evertor muscle weakness represents an important factor affecting chronic ankle instability. For research purposes, ankle evertor strength is assessed by means of isokinetic evaluations. However, this methodology is constraining for daily clinical use. The present study proposes to assess ankle evertor muscle weakness using a new procedure, one that is easily accessible for rehabilitation specialists. To do so, we compared weight bearing ankle inversion control between patients suffering from chronic ankle instability and healthy subjects.

**Methods:** 12 healthy subjects and 11 patients suffering from chronic ankle instability conducted repetitions of one leg weight bearing ankle inversion on a specific ankle destabilization device equipped with a goniometer. Ankle inversion control was performed by means of an eccentric recruitment of evertor muscles. Instructions were to perform, as slow as possible, the ankle inversion while resisting against full body weight applied on the tested ankle.

**Results:** Data clearly showed higher angular inversion velocity peaks in patients suffering from chronic ankle instability. This illustrates an impaired control of weight bearing ankle inversion and, by extension, an eccentric weakness of evertor muscles.

**Interpretation:** The present study supports the hypothesis of a link between the decrease of ankle joint stability and evertor muscle weakness. Moreover, it appears that the new parameter is of use in a clinical setting.

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N. Forestier<sup>1</sup>  
P. Tosch<sup>2</sup>

## The Effects of an Ankle Destabilization Device on Muscular Activity while Walking

Orthopedics & Biomechanics

### Abstract

Chronic Ankle Instability sprain causes are unclear and many factors or mechanisms may contribute to recurrence of this injury. The aim of the study was to investigate how an ankle destabilization device affects the EMG patterns of the ankle muscles during ankle stabilization against inversion. The left foot was equipped with a mechanical device mounted under the heel of the shoe. This mechanical device induces subtalar joint destabilization necessitating the control of ankle muscles. Surface electrodes were placed over the tibialis anterior, the peroneus longus, the peroneus brevis, the gastrocnemius lateral, and the gastrocnemius medial. Nine healthy subjects (mean age  $37 \pm 12$  yr; mean mass  $68 \pm 17$  kg; mean height  $1.73 \pm 0.07$  m) were instructed to walk

normally along a tape fixed on the floor. The ankle destabilization device altered the walking pattern of all subjects. More specifically, the walking pattern is disturbed resulting in higher amplitude of the EMG activity of the peroneal muscles and the Tibialis Anterior and anticipatory reactions in the peroneal muscles. The results suggest that the ankle destabilization device could be beneficial for rehabilitation programs especially during the training of walking. Using this material may help to a specific reinforcement of muscles involved in anti-inversion ankle movement.

### Key words

Electromyography · orthosis

## [ RESEARCH REPORT ]

LUKE DONOVAN, PhD, ATC<sup>1</sup> • JOSEPH M. HART, PhD, ATC<sup>2</sup>  
JAY HERTZEL, PhD, ATC<sup>2</sup>

## Effects of 2 Ankle Destabilization Devices on Electromyography Measures During Functional Exercises in Individuals With Chronic Ankle Instability

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*Therapeutic strategy for external ankle sprains: A preliminary clinical study on the Myolux™ device's efficiency*

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The screenshot shows an IDE with the following components:

- Project Explorer (Left):** Displays the project structure for 'MyoluxMedikEvolution'. It includes packages like 'adapter', 'dao', 'dialog', 'fragment', 'helper', 'model', 'view', and 'res'. The 'dao' package contains classes like 'AbstractDAO', 'AdminDAO', 'CurrentDAO', 'PatientDAO', 'PrescriptionDAO', 'ResultDAO', and 'SeanceDAO'.
- Code Editor (Center):** Shows the implementation of the 'SeanceDAO' class. The visible code includes:
 

```

      addTitleUnderSection(context.getResources().getString(temp[1]));

      List list = new
      if (averages[1]
      list.add(add
      } else if (avera
      list.add(add
      } else if (avera
      list.add(add
      } else {
      list.add(add
      }

      private ListItem addComm
      ListItem item = new
      item.setAlignment(30
      Phrase r = new Phras
      item.add(r);
      r = new Phrase(cont
      item.add(r);
      return item;

      private void addCommentB
      if (exerciseUnder25.
      addTitleUnderSec
      addExerciseUnder

      if (exerciseUnderBet
      addTitleUnderSec
      
```
- Right Pane:** Shows the 'SeanceDAO' class with the following code:
 


```

      ContentValues contentValues = new ContentValues();
      contentValues.put(Patient.COLUMN_NAME_LASTDATE, EditDateView.toFormattedDate(seance.getDate().getTime()));
      PatientDAO.get(context).update(contentValues, Patient._ID + " = " + prescription.getPatientId());
      return seance;

      @Override
      public void delete(int id) {
      int prescriptionId = SeanceDAO.get(context).findOne(id).getPrescriptionId();
      Prescription p = PrescriptionDAO.get(context).findOne(prescriptionId);
      if (p != null) {
      int patientId = p.getPatientId();
      Seance temp = null;
      for (Prescription prescription : PatientDAO.get(context).findOne(patientId).getPrescriptions()) {
      if (!prescription.getSeances().isEmpty()) {
      Seance seance = prescription.getSeances().get(prescription.getSeances().size() - 1);
      if (temp == null) {
      temp = seance;
      } else {
      if (temp.getDate().compareTo(seance.getDate()) < 0) {
      temp = seance;
      }
      }
      }
      }
      ContentValues values = new ContentValues();
      if (temp == null) {
      values.put(Patient.COLUMN_NAME_LASTDATE, "");
      } else {
      values.put(Patient.COLUMN_NAME_LASTDATE, EditDateView.toFormattedDate(temp.getDate().getTime()));
      }
      PatientDAO.get(context).update(values, Patient._ID + " = " + patientId);
      super.delete(id);
      Current current = CurrentDAO.get(context).findOne();
      if (current != null && current.getIdSeance() == id) {
      
```
- Status Bar (Bottom):** Shows 'Gradle build finished in 46s 788ms (25 minutes ago)' and '17:14 CRLF: UTF-8: Git: develop Context: <no context>'.

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
The screenshot shows the home screen of the MYOLUX MEDIK E-VOLUTION application. At the top, there is a grey header bar with a hamburger menu icon on the left, the text "MYOLUX MEDIK E-VOLUTION" in the center, and a "CONNEXION" button with a lock icon on the right, which is circled in orange. Below the header, the main content area features a large hexagonal diagram with six numbered nodes (1-6) and a central orange hexagon. The nodes are labeled as follows:


- Node 1: PROPRIOCEPTION REPOSITIONNEMENT
- Node 2: PROPRIOCEPTION DISSOCIATION
- Node 3: PROPRIOCEPTION FONCTIONNELLE
- Node 4: RENFORCEMENT EXCENTRIQUE
- Node 5: RENFORCEMENT CONCENTRIQUE
- Node 6: RENFORCEMENT FONCTIONNEL

At the bottom right, there is a legend for foot selection:

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System status icons for Bluetooth, Wi-Fi, battery, and silent mode are visible in the top right corner.

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The screenshot displays the 'FICHE PATIENT' (Patient Card) interface. At the top, there is a back arrow, the title 'FICHE PATIENT', and a 'DÉCONNEXION' (Logout) button. Below the title, there is a search bar 'Rechercher...' and a search icon. To the right of the search bar are icons for Bluetooth, a tablet, a battery level indicator, and a mute icon. Below the search bar, there is a 'Trier par' (Sort by) section with radio buttons for 'Nom' (selected) and 'Date'. The patient's name is split into 'Prénom' (Jacky) and 'Nom' (Michel), with a 'Date de naissance' field. A patient card for 'Jacky MICHEL' (10 mars 2016) is shown with a trash icon. Below this is a 'Commentaire' section. A 'Prescription 1' section is highlighted in teal, containing a table of sessions. A button '+ AJOUTER UNE SÉANCE' is circled in orange. The table lists three sessions, each with a refresh icon, a date of '10 mars 2016', and a trash icon. At the bottom, there are two buttons circled in orange: '+ AJOUTER UN PATIENT' on the left and '+ AJOUTER UNE PRÉSCRIPTION' on the right. The bottom of the screen shows the Android navigation bar.

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The screenshot displays the MYOLUX MEDIK E-EVOLUTION application interface. At the top, there is a navigation bar with a hamburger menu icon on the left, the text "MYOLUX MEDIK E-EVOLUTION" in the center, and a "CONNEXION" button with a lock icon on the right. Below the navigation bar, the title "PROPRIOCEPTION FONCTIONNELLE" is centered. The main content area features a hexagonal diagram with six numbered nodes (1-6) and a central orange hexagon. The nodes are labeled as follows:

- 1: PROPRIOCEPTION REPOSITIONNEMENT
- 2: PROPRIOCEPTION DISSOCIATION
- 3: (unlabeled)
- 4: RENFORCEMENT EXCENTRIQUE
- 5: RENFORCEMENT CONCENTRIQUE
- 6: RENFORCEMENT FONCTIONNEL

At the bottom right, there is a legend for foot selection:

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System status icons for Bluetooth, Wi-Fi, mobile data, battery, and volume are visible in the top right corner.

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
The screenshot shows the app interface for a proprioception test. At the top, there is a navigation bar with a back arrow, the title "PROPRIOCEPTION REPOSITIONNEMENT", and buttons for "FICHE PATIENT" and "DÉCONNEXION". Below the navigation bar, there are icons for Bluetooth, a tablet, a battery level indicator, and a mute icon. On the left side, there is an image of a foot wearing a blue and white sensor device. In the center, there is a button labeled "REVOIR LE DÉTAIL DES INSTRUCTIONS". Below this, there is a text instruction: "Assis, avec l'avant-pied décollé du sol et l'arrière pied sur un support d'environ 2.5 cm." There are three selection options: "Choisir le pied" with radio buttons for "Pied Gauche" (selected) and "Pied Droit"; "Choisir la position" with radio buttons for "Assis" (selected), "Sur deux pieds", and "Sur un pied"; and "Choisir la difficulté" with radio buttons for "Facile", "Normal", and "Difficile". At the bottom, there is a horizontal slider with a scale from -30 to 5. A green circle is positioned at approximately -25, and an orange circle is positioned at approximately 5. On the left side of the slider, there is a green circular button labeled "Arrêter le TEST" with a white square icon and a white circle containing the number "07".

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The screenshot shows the application interface for the 'PROPRIOCEPTION REPOSITIONNEMENT' test. At the top, there is a back arrow, an information icon, and the title 'PROPRIOCEPTION REPOSITIONNEMENT'. To the right are buttons for 'FICHE PATIENT' and 'DÉCONNEXION'. Below the title, there are icons for Bluetooth, screen rotation, battery, and volume. A central button reads 'REVOIR LE DÉTAIL DES INSTRUCTIONS'. An image of a foot on a MyoLux device is shown on the left. The instructions state: 'Assis, avec l'avant-pied décollé du sol et l'arrière pied sur un support d'environ 2.5 cm.' Below this are three selection options: 'Choisir le pied : Pied Gauche (selected) / Pied Droit', 'Choisir la position : Assis (selected) / Sur deux pieds / Sur un pied', and 'Choisir la difficulté : Facile / Normal / Difficile'. A large green circular button says 'Lancer le TEST' with a play icon and a '10' timer. A horizontal slider at the bottom ranges from -30 to 5, with an orange knob currently positioned at -15. The bottom of the screen shows standard Android navigation icons.

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© 2014 Elsevier Ltd. All rights reserved.